


APPLICATION NOTE

FMBIO® Detection of SYPRO® Red Dye Stained Proteins in Polyacrylamide Gels



1 **Two-fold dilution series of human serum stained with SYPRO® Red dye.** Serum was serially diluted (Lane 1 = 13.0 µg) and loaded onto a 4-20%, 1-mm Tris-Glycine gel (NOVEX) and electrophoresed for 90 minutes in a 1X Tris Glycine running buffer¹ with 0.05% SDS. Gels were stained in SYPRO® Red dye (1:5000 in 7.5% acetic acid) with gentle agitation. To reduce background, a destaining step was performed in which gels were gently agitated in 7.5% acetic acid for 45 minutes. The image was acquired on the FMBIO® II Fluorescence Imaging System using a 625/15 nm filter at 100% photomultiplier tube sensitivity, with 256 repeats and a spatial resolution of 300 x 300 dots per inch.

Applications

SYPRO® Red protein gel stain was developed to provide fast, simple, sensitive staining of proteins in electrophoretic gels. The FMBIO® II Fluorescence Imaging System can be used to image the following types of protein gels stained with SYPRO Red dye:

- SDS-denaturing Gels
- Native Gels
- 2-D Gels

Benefits


Scanning SYPRO Red dye stained gels on the FMBIO II Fluorescence Imaging System provides a number of benefits compared to Coomassie® Brilliant Blue and silver staining^{2,3}.

- Staining does not modify antigenic sites, so they are left available for antibody binding.
 - Polypeptide sequence and degree of glycosylation do not affect staining intensity and quality.
 - Detection sensitivities are comparable to silver staining and are higher than Coomassie staining.
 - Nucleic acids found in protein mixtures from cell or tissue extracts are not stained.
 - SYPRO Red stain is optimally matched to the FMBIO laser for maximum imaging sensitivity, with minimum background fluorescence.
 - The combination of the FMBIO and SYPRO Red stain results in greater dynamic range than instant film or CCD cameras.
 - FMBIO Analysis Software allows a wide variety of quantitation capabilities.
 - Digital images can be created, published, and archived, eliminating the need to dry and save gels.
- Data acquisition with the FMBIO is fast and easy.
 - Gel fixation is unnecessary, resulting in efficient transfer of protein bands to membranes for western blot analysis.

Sensitivity
625 pg/band

Dynamic Range
2.2 orders of magnitude (625 pg/band to 100 ng/band)

Dye Spectra



Tips for FMBIO Imaging

1. Many factors contribute to the success of protein detection. Variables include: gel type and concentration, sample preparation, band size, dye concentration, staining time, destaining time, gray level settings, and laser focusing point. Further optimization of protocols may result in higher detection sensitivity.
2. For best scanning results, avoid dust specks by using only powderless gloves, rinsing gloves with distilled water, thoroughly cleaning all containers used for staining, and filtering all buffers and solutions with a 0.45µm filter.
3. For lowest background, place the gel on a low-fluorescence glass plate or sample holder when scanning. Do not scan on plastic wrap unless the wrap has been tested for low background fluorescence.
4. Keep solutions containing SYPRO Red stain covered by aluminum foil to prevent photobleaching of the dye

by ambient light. Staining reagent diluted in buffer or acetic acid and stored protected from light at 4°C may be reused for at least three months.

5. Before opening each vial, warm it to room temperature and then briefly centrifuge. If dye particles are present, they should be redissolved by briefly vortexing the tube after warming.

6. For detailed protocols on the usage, handling and disposal of SYPRO Red stain, see Molecular Probes Product Information Sheet TD005⁴.

FMBIO Analysis Software provides easy-to-use features for image analysis:

- Band Sizing
- 1D Volume Calculation
- 1D Peak Height
- Lane Traces
- Spot Quantitation
- Background Subtraction
- Standard Curve
- Text and Other Annotation
- Gray Level Adjustment
- User-defined Scan Options

References

1. Ausubel, F.M. *et al.* (1992). Short Protocols in Molecular Biology, 2nd Ed., John Wiley and Sons, New York.
2. SYPRO® Protein Gel Stains. (1996). Molecular Probes Product Information Sheet MP 6650.
3. Ishino Y., *et al.* (1992). Practical Applications in Molecular Biology of Sensitive Fluorescence Detection by a Laser-excited Fluorescence Imager Analyzer. *BioTechnique* **13**, 936.
4. SYPRO® Protein Gel Stains (1997). Molecular Probes Product Information Sheet TD 005.



Hitachi Software
Engineering America, Ltd.

601 Gateway Blvd., Suite 500
South San Francisco, CA 94080
Tel: 800-624-6176 or 650-615-9600
Fax: 650-615-7699

Internet: <http://www.hitachi-soft.com/gs>
Email: gene@hitsoft.com

Hitachi Software
Engineering Co., Ltd.

5-79 Onoe-Cho, Naka-ku
Yokohama 231-0015, Japan
Tel: 81-45-681-2252
Fax: 81-45-681-4399

Hitachi Software
Engineering Europe S.A.

Parc de Limère-Z.I. d'Ardon-B.P. 629
45166 OLIVET Cedex-France
Tel: (+33) 238-69-86-90
Fax: (+33) 238-69-86-99

Email: gene@software.hitachi.fr

FMBIO is a registered trademark of Hitachi Software Engineering Co., Ltd. SYPRO is a registered trademark of Molecular Probes, Inc. NOVEX is a trademark of Novel Experimental Technology. Coomassie is a registered trademark of Imperial Chemical Industries, Ltd.

©1998 Hitachi Software Engineering America, Ltd. All rights reserved.