

SAMPLE PROTOCOL FOR WASHED CAPTURE SANDWICH IMMUNOASSAY USING MAGNETIC MICROSPHERES

Microspheres should be protected from prolonged exposure to light throughout this procedure.

1. Select the appropriate antibody-coupled microsphere sets.
2. Resuspend the microspheres by vortex and sonication for approximately 20 seconds.
3. Prepare a Working Microsphere Mixture by diluting the coupled microsphere stocks to a final concentration of 100 microspheres of each set/ μL in PBS-1% BSA. (Note: 50 μL of Working Microsphere Mixture is required for each reaction.) See **Technical Note 1**.
4. Aliquot 50 μL of the Working Microsphere Mixture into the appropriate wells of a round-bottom plate.
5. Add 50 μL of PBS-1% BSA to each background well.
6. Add 50 μL of standard or sample to the appropriate wells.
7. Mix the reactions gently by pipetting up and down several times with a multi-channel pipettor.
8. Cover the plate and incubate for 30 minutes at room temperature on a plate shaker set to approximately 800 rpm.
9. Place the plate into the magnetic separator and allow separation to occur for 30 to 60 seconds. See **Technical Note 2**.
10. Use a multi-channel pipette to carefully aspirate the supernatant from each well. Take care not to disturb the microspheres.
11. Leave the plate in the magnetic separator for the following wash steps:
 - a. Add 100 μL PBS-1% BSA to each well.
 - b. Use a multi-channel pipette to carefully aspirate the supernatant from each well. Take care not to disturb the microspheres.
 - c. Repeat steps **a.** and **b.** above.
12. Remove the plate from the magnetic separator and resuspend the microspheres in 50 μL of PBS-1% BSA by gently pipetting up and down several times using a multi-channel pipettor.
13. Dilute the biotinylated detection antibody to 4 $\mu\text{g}/\text{mL}$ in PBS-1% BSA. (Note: 50 μL of diluted detection antibody is required for each reaction.) See **Technical Note 3**.
14. Add 50 μL of the diluted detection antibody to each well.
15. Mix the reactions gently by pipetting up and down several times with a multi-channel pipettor.
16. Cover the plate and incubate for 30 minutes at room temperature on a plate shaker set to approximately 800 rpm.
17. Place the plate into the magnetic separator and allow separation to occur for 30 to 60 seconds.

18. Use a multi-channel pipette to carefully aspirate the supernatant from each well. Take care not to disturb the microspheres.
19. Leave the plate in the magnetic separator for the following wash steps:
 - a. Add 100 μ L PBS-1% BSA to each well.
 - b. Use a multi-channel pipette to carefully aspirate the supernatant from each well. Take care not to disturb the microspheres.
 - c. Repeat steps **a.** and **b.** above.
20. Remove the plate from the magnetic separator and resuspend the microspheres in 50 μ L of PBS-1% BSA by gently pipetting up and down several times with a multi-channel pipettor.
21. Dilute streptavidin-R-phycoerythrin reporter to 4 μ g/mL in PBS-1% BSA. (Note: 50 μ L of diluted streptavidin-R-phycoerythrin is required for each reaction.) See **Technical Note 3**.
22. Add 50 μ L of the diluted streptavidin-R-phycoerythrin to each well.
23. Mix the reactions gently by pipetting up and down several times with a multi-channel pipettor.
24. Cover the plate and incubate for 30 minutes at room temperature on a plate shaker set to approximately 800 rpm.
25. Place the plate into the magnetic separator and allow separation to occur for 30 to 60 seconds.
26. Use a multi-channel pipette to carefully aspirate the supernatant from each well. Take care not to disturb the microspheres.
27. Leave the plate in the magnetic separator for the following wash steps:
 - a. Add 100 μ L PBS-1% BSA to each well.
 - b. Use a multi-channel pipette to carefully aspirate the supernatant from each well. Take care not to disturb the microspheres.
 - c. Repeat steps **a.** and **b.** above.
28. Remove the plate from the magnetic separator and resuspend the microspheres in 100 μ L of PBS-1% BSA by gently pipetting up and down several times with a multi-channel pipettor.
29. Analyze 50-75 μ L on the Luminex analyzer according to the system manual.

Technical Note 1: Either PBS-1% BSA or PBS-BN (PBS, 1% BSA, 0.05% Azide, pH 7.4) may be used as Assay Buffer.

Technical Note 2: For a list of magnetic separator plates, see **Recommended Materials for Magnetic Microspheres**. Optimal separation time may vary with the type of separator used.

Technical Note 3: Concentrations should be optimized for specific reagents, assay conditions, level of multiplexing, etc. in use.